

IFNAR1 Mouse mAb

Catalog No: #63766

Package Size: #63766-1 50ul #63766-2 100ul


Description

Product Name	IFNAR1 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;IHC;FC
Species Reactivity	Human,Mouse
Immunogen Description	Purified recombinant fragment of human IFNAR1 (AA: 28-156) expressed in E. Coli.
Target Name	IFNAR1
Other Names	AVP; IFRC; IFNAR; IFNBR; IFN-alpha-REC
Accession No.	P17181
Calculated MW	63.5kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000IHC:1/200 - 1/1000FC:1/200 - 1/400ELISA:1/10000


Images




Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)




Western blot analysis using IFNAR1 mouse mAb against mouse lung (1) lysate.



Flow cytometric analysis of K562 cells using IFNAR1 mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human kidney cancer tissues using IFNAR1 mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human rectum cancer tissues using IFNAR1 mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.