

IGF1R Mouse mAb

Catalog No: #63768

Package Size: #63768-1 50ul #63768-2 100ul


Description

Product Name	IGF1R Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG2a
Applications	WB;FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human IGF1R (AA: extra(741-935)) expressed in E. Coli.
Target Name	IGF1R
Other Names	IGFR; CD221; IGFIK; JTK13
Accession No.	P08069
Calculated MW	154.8kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000FC:1/200 - 1/400ELISA:1/10000


Images




Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Western blot analysis using IGF1R mouse mAb against NIH3T3 (1), cell lysate.



Flow cytometric analysis of HepG2 cells using IGF1R mouse mAb (green) and negative control (red).



Flow cytometric analysis of Lovo cells using IGF1R mouse mAb (green) and negative control (red).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.