

IGFBP2 Mouse mAb

Catalog No: #63771

Package Size: #63771-1 50ul #63771-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	IGFBP2 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	IHC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of protein IGFBP2 (aa180-328) expressed in E. Coli.
Target Name	IGFBP2
Other Names	IBP2; IGF-BP53
Accession No.	P18065
Formulation	Purified antibody in PBS with 0.05% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

IHC:1/200 - 1/1000ELISA:1/10000

Images



Immunohistochemical analysis of paraffin-embedded human lung cancer (A), rectum (B), prostate (C), colon cancer (D) showing cytoplasmic localization using IGFBP2 mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.