

IL10RB Mouse mAb

Catalog No: #63776

Package Size: #63776-1 50ul #63776-2 100ul



Description


Product Name	IL10RB Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;FC
Species Reactivity	Human,Rat
Immunogen Description	Purified recombinant fragment of human IL10RB (AA: extra 20-220) expressed in E. Coli.
Target Name	IL10RB
Other Names	CRFB4; CRF2-4; D21S58; D21S66; CDW210B; IL-10R2
Accession No.	Q08334
Calculated MW	37kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

WB:1/500 - 1/2000FC:1/200 - 1/400ELISA:1/10000

Images





Flow cytometric analysis of HL-60 cells using IL10RB mouse mAb (green) and negative control (red).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.