

IL1R1 Mouse mAb

Catalog No: #63780

Package Size: #63780-1 50ul #63780-2 100ul


Description

Product Name	IL1R1 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human IL1R1 (AA: 18-167) expressed in E. Coli.
Target Name	IL1R1
Other Names	P80; IL1R; IL1RA; CD121A; D2S1473; IL-1R-alpha
Accession No.	P14778
Calculated MW	65.4kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000FC:1/200 - 1/400ELISA:1/10000


Images



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Western blot analysis using IL1R1 mouse mAb against Hela (1) cell lysate.



Flow cytometric analysis of HeLa cells using IL1R1 mouse mAb (green) and negative control (red).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.