

IL28A Mouse mAb

Catalog No: #63783

Package Size: #63783-1 50ul #63783-2 100ul


Description

|                       |                                                                                          |
|-----------------------|------------------------------------------------------------------------------------------|
| Product Name          | IL28A Mouse mAb                                                                          |
| Host Species          | Mouse                                                                                    |
| Clonality             | Monoclonal                                                                               |
| Isotype               | Mouse IgG1                                                                               |
| Applications          | WB;IHC                                                                                   |
| Species Reactivity    | Human                                                                                    |
| Immunogen Description | Purified recombinant fragment of human IL28A (AA: 1-200) expressed in E. Coli.           |
| Target Name           | IL28A                                                                                    |
| Other Names           | IFNL2; IL-28A                                                                            |
| Accession No.         | Q8IZJ0                                                                                   |
| Calculated MW         | 22.3kDa                                                                                  |
| Formulation           | Purified antibody in PBS with 0.05% sodium azide                                         |
| Storage               | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles. |


Application Details

WB:1/500 - 1/2000IHC:1/200 - 1/1000ELISA:1/10000


Images




Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Western blot analysis using IL28A mouse mAb against SK-Br-3 (1) cell lysate.



Immunohistochemical analysis of paraffin-embedded human brain tissues using IL28A mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human esophageal cancer tissues using IL28A mouse mAb with DAB staining.

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.