

ITGA5 Mouse mAb

Catalog No: #63796

Package Size: #63796-1 50ul #63796-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	ITGA5 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG2a
Applications	IHC;FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human ITGA5 expressed in E. Coli.
Target Name	ITGA5
Other Names	FNRA; CD49e; VLA5A; ITGA5
Accession No.	P08648
Calculated MW	114kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

IHC:1/200 - 1/1000FC:1/200 - 1/400ELISA:1/10000

Images



Flow cytometric analysis of Hela cells using ITGA5 mouse mAb (green) and negative control (purple).



Immunohistochemical analysis of paraffin-embedded human Placenta tissues using ITGA5 mouse mAb

Note: This product is for in vitro research use only and is not intended for use in humans or animals.