

ITGB4 Mouse mAb

Catalog No: #63798

Package Size: #63798-1 50ul #63798-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	ITGB4 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	IHC;FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human ITGB4 (AA: 1619-1822) expressed in E. Coli.
Target Name	ITGB4
Other Names	CD104
Accession No.	P16144
Calculated MW	202kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

IHC:1/200 - 1/1000FC:1/200 - 1/400ELISA:1/10000


Images




Red: Control Antigen (100ng); Purple:Antigen (10ng); Green: Antigen (50ng);Blue: Antigen (100ng);



Flow cytometric analysis of A549 cells using ITGB4 mouse mAb (green) and negative control (purple).



Immunohistochemical analysis of paraffin-embedded cervical cancer tissues using ITGB4 mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded esophageal cancer tissues using ITGB4 mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.