

JAM3 Mouse mAb

Catalog No: #63803

Package Size: #63803-1 50ul #63803-2 100ul



SAB[®]
Signalway Antibody

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	JAM3 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human JAM3 (AA: extra 32-241) expressed in E. Coli.
Target Name	JAM3
Other Names	JAMC; JAM-2; JAM-3; JAM-C
Accession No.	Q9BX67
Calculated MW	35kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000FC:1/200 - 1/400ELISA:1/10000


Images



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Western blot analysis using JAM3 mouse mAb against SH-SY5Y (1) cell lysate.



Flow cytometric analysis of Hek293 cells using JAM3 mouse mAb (green) and negative control (red).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.