

JUP Mouse mAb

Catalog No: #63804

Package Size: #63804-1 50ul #63804-2 100ul


Description

Product Name	JUP Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;IHC;IF
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human JUP expressed in E. Coli.
Target Name	JUP
Other Names	DP3; PDGB; PKGB; CTNNG; DPIII; ARVD12
Accession No.	P14923
Calculated MW	82kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000IHC:1/200 - 1/1000ICC:1/200 - 1/1000ELISA:1/10000


Images




Red: Control Antigen (100ng); Purple: Antigen (10ng); Green: Antigen (50ng); Blue: Antigen (100ng);




Western blot analysis using JUP mouse mAb against T47D (1), MCF-7 (2), SKBR-3 (3), A431 (4) and HEK293 (5) cell lysate.



Immunofluorescence analysis of U251 cells using JUP mouse mAb (green). Blue: DRAQ5 fluorescent DNA dye.



Immunohistochemical analysis of paraffin-embedded human rectum tissues using JUP mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human stomach cancer tissues using JUP mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.