

KARS Mouse mAb

Catalog No: #63805

Package Size: #63805-1 50ul #63805-2 100ul


Description

Product Name	KARS Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	IHC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of KARS(aa90-174) expressed in E. Coli.
Target Name	KARS
Other Names	KRS; KARS2; KIAA0070
Accession No.	Q15046
Formulation	Ascitic fluid containing 0.03% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

IHC:1/200 - 1/1000ELISA:1/10000

Images



Immunohistochemical analysis of paraffin-embedded human cervical carcinoma, showing cytoplasmic localization with DAB staining using KARS mouse mAb.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.