

KBTBD8 Mouse mAb

Catalog No: #63808

Package Size: #63808-1 50ul #63808-2 100ul


Description

Product Name	KBTBD8 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human KBTBD8 (AA: 264-464) expressed in E. Coli.
Target Name	KBTBD8
Other Names	TA-KRP
Accession No.	Q8NFY9
Calculated MW	68.8kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000FC:1/200 - 1/400ELISA:1/10000


Images



Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);



Western blot analysis using KBTBD8 mouse mAb against A431 (1), Jurkat (2), Hela (3), K562 (4), and HEK293 (5) cell lysate.



Flow cytometric analysis of A431 cells using KBTBD8 mouse mAb (green) and negative control (red).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.