

KCNK1 Mouse mAb

Catalog No: #63810

Package Size: #63810-1 50ul #63810-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	KCNK1 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG2b
Applications	IHC;FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human KCNK1 (AA: extra mix) expressed in E. Coli.
Target Name	KCNK1
Other Names	DPK; HOHO; K2P1; KCNO1; TWIK1; K2p1.1; TWIK-1
Accession No.	O00180
Calculated MW	38.1kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

IHC:1/200 - 1/1000FC:1/200 - 1/400ELISA:1/10000


Images




Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)




Flow cytometric analysis of C6 cells using KCNK1 mouse mAb (green) and negative control (red).



Flow cytometric analysis of HepG2 cells using KCNK1 mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human esophageal cancer tissues using KCNK1 mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human stomach cancer tissues using KCNK1 mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.