

KDM1A Mouse mAb

Catalog No: #63813

Package Size: #63813-1 50ul #63813-2 100ul


Description

Product Name	KDM1A Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;IHC;IF;FC
Species Reactivity	Human,Mouse,Monkey
Immunogen Description	Purified recombinant fragment of human KDM1A (AA: 709-876) expressed in E. Coli.
Target Name	KDM1A
Other Names	AOF2; CPRF; KDM1; LSD1; BHC110
Accession No.	O60341
Calculated MW	93kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000IHC:1/200 - 1/1000ICC:1/500 - 1/2000FC:1/200 - 1/400ELISA:1/10000


Images




Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)




Western blot analysis using KDM1A mouse mAb against COS7 (1), NIH/3T3 (2), and A549 (3) cell lysate.



Flow cytometric analysis of Hela cells using KDM1A mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human cervical cancer tissues using KDM1A mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human esophageal cancer tissues using KDM1A mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.