

KDM3A Mouse mAb

Catalog No: #63814

Package Size: #63814-1 50ul #63814-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	KDM3A Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;IHC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human KDM3A expressed in E. Coli.
Target Name	KDM3A
Other Names	TSGA; JMJD1; JHDM2A; JHMD2A; JMJD1A; KIAA0742; DKFZp686A24246; DKFZp686P07111; KDM3A
Accession No.	Q9Y4C1
Calculated MW	147kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000IHC:1/200 - 1/1000ELISA:1/10000

Images



Western blot analysis using KDM3A mouse mAb against HeLa (1) and HepG2 (2) cell lysate.



Immunohistochemical analysis of paraffin-embedded human colonic cancer tissues (left) and lung cancer tissues (right) using KDM3A mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.