

KDM4A Mouse mAb

Catalog No: #63817

Package Size: #63817-1 50ul #63817-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	KDM4A Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	IF;FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human KDM4A (AA: 932-1057) expressed in E. Coli.
Target Name	KDM4A
Other Names	JMJD2; JHDM3A; JMJD2A; TDRD14A
Accession No.	O75164
Calculated MW	120.6kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

ICC:1/200 - 1/1000FC:1/200 - 1/400ELISA:1/10000


Images



Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);



Immunofluorescence analysis of HeLa cells using KDM4A mouse mAb. Blue: DRAQ5 fluorescent DNA dye. Red: Actin filaments have been labeled with Alexa Fluor- 555 phalloidin.



Immunofluorescence analysis of HeLa cells using KDM4A mouse mAb (green). Blue: DRAQ5 fluorescent DNA dye. Red: Actin filaments have been labeled with Alexa Fluor- 555 phalloidin. Secondary antibody from Fisher (Cat#: 35503)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.