

KDM4B Mouse mAb

Catalog No: #63818

Package Size: #63818-1 50ul #63818-2 100ul


Description

Product Name	KDM4B Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	FC
Species Reactivity	Human,Mouse
Immunogen Description	Purified recombinant fragment of human KDM4B (AA: 960-1096) expressed in E. Coli.
Target Name	KDM4B
Other Names	JMJD2B; TDRD14B
Accession No.	O94953
Calculated MW	122kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

FC:1/200 - 1/400ELISA:1/10000

Images



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Flow cytometric analysis of NIH/3T3 cells using KDM4B mouse mAb (green) and negative control (red).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.