

KEAP1 Mouse mAb

Catalog No: #63820

Package Size: #63820-1 50ul #63820-2 100ul


Description

Product Name	KEAP1 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;IF;FC
Species Reactivity	Human,Mouse
Immunogen Description	Purified recombinant fragment of human KEAP1 (AA: 380-624) expressed in E. Coli.
Target Name	KEAP1
Other Names	INrf2; KLHL19
Accession No.	Q14145
Calculated MW	69.7kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

ICC:1/50FC:1/200 - 1/400ELISA:1/10000


Images




Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);



Western blot analysis using KEAP1 mouse mAb against NIH3T3 (1), and A549 (2) cell lysate.



Immunofluorescence analysis of HeLa cells using KEAP1 mouse mAb (green). Blue: DRAQ5 fluorescent DNA dye. Red: Actin filaments have been labeled with Alexa Fluor-555 phalloidin.



Flow cytometric analysis of HepG2 cells using KEAP1 mouse mAb (green) and negative control (purple).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.