

KIR3DL1 Mouse mAb

Catalog No: #63822

Package Size: #63822-1 50ul #63822-2 100ul



Description


Product Name	KIR3DL1 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG2b
Applications	WB;FC
Species Reactivity	Human,Rat
Immunogen Description	Purified recombinant fragment of human KIR3DL1 (AA: extra 22-340) expressed in HEK293 cells.
Target Name	KIR3DL1
Other Names	KIR; NKB1; NKAT3; NKB1B; NKAT-3; CD158E1; KIR3DL1/S1
Accession No.	P43629
Calculated MW	49kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

WB:1/500 - 1/2000FC:1/200 - 1/400ELISA:1/10000

Images





Flow cytometric analysis of Jurkat cells using KIR3DL1 mouse mAb (green) and negative control (red).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.