

KIT Mouse mAb

Catalog No: #63824

Package Size: #63824-1 50ul #63824-2 100ul


Description

Product Name	KIT Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;IHC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human KIT expressed in E. Coli.
Target Name	KIT
Other Names	PBT; SCFR; C-Kit; CD117
Accession No.	P10721
Calculated MW	145kDa
Formulation	Ascitic fluid containing 0.03% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000IHC:1/200 - 1/1000ELISA:1/10000


Images



Red: Control Antigen (100ng); Purple: Antigen (10ng); Green: Antigen (50ng); Blue: Antigen (100ng);



Western blot analysis using KIT mouse mAb against Jurkat (1) and Hela (2) cell lysate.



Immunohistochemical analysis of paraffin-embedded human gastrointestinal stromal tumor using KIT mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.