

KLF15 Mouse mAb

Catalog No: #63827

Package Size: #63827-1 50ul #63827-2 100ul


Description

Product Name	KLF15 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human KLF15 expressed in E. Coli.
Target Name	KLF15
Other Names	KKLF; DKFZp779M1320
Accession No.	Q9UIH9
Calculated MW	44kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

WB:1/500 - 1/2000ELISA:1/10000

Images



Western blot analysis using KLF15 mouse mAb against HepG2 (1) and SMMC-7721 (2) cell lysate.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.