

KLF4 Mouse mAb

Catalog No: #63829

Package Size: #63829-1 50ul #63829-2 100ul


Description

Product Name	KLF4 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	IHC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human KLF4 (AA: 2-470) expressed in E. Coli.
Target Name	KLF4
Other Names	EZF; GKLF; KLF4
Accession No.	O43474
Calculated MW	55kDa
Formulation	Ascitic fluid containing 0.03% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

IHC: 1/200 - 1/1000ELISA: 1/10000

Images



Immunohistochemical analysis of paraffin-embedded human lung cancer tissues (left) and human rectum tissues (right) using KLF4 mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.