

KLHL11 Mouse mAb

Catalog No: #63832

Package Size: #63832-1 50ul #63832-2 100ul


Description

Product Name	KLHL11 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;IHC;IF
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human KLHL11 expressed in E. Coli.
Target Name	KLHL11
Other Names	FLJ10572; KLHL11
Accession No.	Q9NVR0
Calculated MW	80kDa
Formulation	Ascitic fluid containing 0.03% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

IHC:1/200 - 1/1000ICC:1/200 - 1/1000ELISA:1/10000


Images



Western blot analysis using KLHL11 mouse mAb against HeLa (1) and MCF-7 (2) cell lysate.



Immunofluorescence analysis of NIH/3T3 cells using KLHL11 mouse mAb (green). Blue: DRAQ5 fluorescent DNA dye. Red: Actin filaments have been labeled with Alexa Fluor-555 phalloidin.



Immunohistochemical analysis of paraffin-embedded human liver cancer (right) and colon cancer tissues (left) using KLHL11 mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.