

KLHL12 Mouse mAb

Catalog No: #63833

Package Size: #63833-1 50ul #63833-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	KLHL12 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human KLHL12 expressed in E. Coli.
Target Name	KLHL12
Other Names	DKIR; C3IP1; FLJ27152; KLHL12
Accession No.	Q53G59
Calculated MW	63kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

WB: 1/500 - 1/2000 ELISA: 1/10000

Images



Western blot analysis using KLHL12 mouse mAb against Hela (1) cell lysate.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.