

KLHL22 Mouse mAb

Catalog No: #63834

Package Size: #63834-1 50ul #63834-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	KLHL22 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;IF
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human KLHL22 expressed in E. Coli.
Target Name	KLHL22
Other Names	KELCHL; KLHL22
Accession No.	Q53GT1
Calculated MW	72kDa
Formulation	Ascitic fluid containing 0.03% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000 ICC:1/200 - 1/1000 ELISA:1/10000

Images



Western blot analysis using KLHL22 mouse mAb against mouse brain tissues lysate.



Immunofluorescence analysis of U251 cells using KLHL22 mouse mAb (green). Blue: DRAQ5 fluorescent DNA dye. Red: Actin filaments have been labeled with Alexa Fluor-555 phalloidin.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.