

KMT2C Mouse mAb

Catalog No: #63838

Package Size: #63838-1 50ul #63838-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	KMT2C Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	IHC;IF;FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human KMT2C (AA: 1-205) expressed in E. Coli.
Target Name	KMT2C
Other Names	HALR; MLL3; KLEFS2
Accession No.	Q8NEZ4
Calculated MW	541.4kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

IHC:1/200 - 1/1000 ICC:1/200 - 1/1000 FC:1/200 - 1/400 ELISA:1/10000


Images




Black line: Control Antigen (100 ng); Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng)



Immunofluorescence analysis of HeLa cells using KMT2C mouse mAb (green). Blue: DRAQ5 fluorescent DNA dye. Red: Actin filaments have been labeled with Alexa Fluor- 555 phalloidin. Secondary antibody from Fisher (Cat#: 35503)



Flow cytometric analysis of HeLa cells using KMT2C mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human cervical cancer tissues using KMT2C mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.