

KMT5A Mouse mAb

Catalog No: #63841

Package Size: #63841-1 50ul #63841-2 100ul


Description

Product Name	KMT5A Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	IHC;FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human KMT5A (AA: 109-272) expressed in E. Coli.
Target Name	KMT5A
Other Names	SET8; SET07; SETD8; PR-Set7; PR/SET07
Accession No.	Q9NQR1
Calculated MW	42.9kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

IHC:1/200 - 1/1000FC:1/200 - 1/400ELISA:1/10000


Images



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Flow cytometric analysis of HeLa cells using KMT5A mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human bladder cancer tissues using KMT5A mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.