

KRT15 Mouse mAb

Catalog No: #63847

Package Size: #63847-1 50ul #63847-2 100ul


Description

Product Name	KRT15 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG2a
Applications	WB;IHC;IF;FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of KRT15 expressed in E. Coli.
Target Name	KRT15
Other Names	K15; CK15; K1CO; KRT1
Accession No.	P19012
Calculated MW	49kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

IHC:1/200 - 1/1000ICC:1/200 - 1/1000FC:1/200 - 1/400ELISA:1/10000


Images




Western blot analysis using KRT15 mouse mAb against A431 cell lysate.



Immunofluorescence analysis of HepG2(left) and PACN-1 (right) cells using KRT15 mouse mAb (green). Red: Actin filaments have been labeled with DY-554 phalloidin. Blue: DRAQ5 fluorescent DNA dye.



Flow cytometric analysis of PACN-1 cells using KRT15 mouse mAb (green) and negative control (purple).



Immunohistochemical analysis of paraffin-embedded human Tonsil tissues using KRT15 mouse mAb

Note: This product is for in vitro research use only and is not intended for use in humans or animals.