

KRT2 Mouse mAb

Catalog No: #63848

Package Size: #63848-1 50ul #63848-2 100ul


Description

Product Name	KRT2 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG2a
Applications	WB;IHC;IF;FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human KRT2 (AA: 177-490) expressed in E. Coli.
Target Name	KRT2
Other Names	K2e; KRTE; CK-2e; KRT2A; KRT2E
Accession No.	P35908
Calculated MW	65.4kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000IHC:1/200 - 1/1000ICC:1/200 - 1/1000FC:1/200 - 1/400ELISA:1/10000


Images




Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)




Western blot analysis using KRT2 mouse mAb against A431 (1) cell lysate.




Immunofluorescence analysis of HeLa cells using KRT2 mouse mAb (green). Blue: DRAQ5 fluorescent DNA dye. Red: Actin filaments have been labeled with Alexa Fluor- 555 phalloidin. Secondary antibody from Fisher (Cat#: 35503)



Flow cytometric analysis of HeLa cells using KRT2 mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human lung cancer tissues using KRT2 mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human bladder cancer tissues using KRT2 mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.