

KRT6A Mouse mAb

Catalog No: #63849

Package Size: #63849-1 50ul #63849-2 100ul


Description

Product Name	KRT6A Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG2b
Applications	WB;IF;FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human KRT6A (AA: 253-352) expressed in E. Coli.
Target Name	KRT6A
Other Names	KRT6A; K6A; K6C; K6D; PC3; CK6C; CK6D; CK-6C; CK-6E; KRT6C; KRT6D
Accession No.	P02538
Calculated MW	60kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000ICC:1/200 - 1/1000FC:1/200 - 1/400ELISA:1/10000


Images




Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Western blot analysis using KRT6A mouse mAb against Hela (1), MCF-7 (2), SW480 (3), HepG2 (4), and T47D (5) cell lysate.



Immunofluorescence analysis of HeLa cells using KRT6A mouse mAb (green). Blue: DRAQ5 fluorescent DNA dye. Red: Actin filaments have been labeled with Alexa Fluor- 555 phalloidin. Secondary antibody from Fisher (Cat#: 35503)



Flow cytometric analysis of HeLa cells using KRT6A mouse mAb (green) and negative control (red).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.