

KRT9 Mouse mAb

Catalog No: #63851

Package Size: #63851-1 50ul #63851-2 100ul


Description

Product Name	KRT9 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG2b
Applications	WB;IHC;FC
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Purified recombinant fragment of human KRT9 (AA: 153-460) expressed in E. Coli.
Target Name	KRT9
Other Names	K9; CK-9; EPPK
Accession No.	P35527
Calculated MW	62kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000IHC:1/200 - 1/1000FC:1/200 - 1/400ELISA:1/10000


Images




Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)




Western blot analysis using KRT9 mouse mAb against HepG2 (1), MCF-7 (2), NIH/3T3 (3), rat heart (4), mouse heart (5), Hela (6), and HUVEC (7) cell lysate.



Flow cytometric analysis of A431 cells using KRT9 mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human bladder cancer tissues using KRT9 mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human esophageal cancer tissues using KRT9 mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.