

KRT9 Mouse mAb

Catalog No: #63852

Package Size: #63852-1 50ul #63852-2 100ul


Description

Product Name	KRT9 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG2b
Applications	WB;FC
Species Reactivity	Human,Monkey
Immunogen Description	Purified recombinant fragment of human KRT9 (AA: 153-460) expressed in E. Coli.
Target Name	KRT9
Other Names	K9; CK-9; EPPK
Accession No.	P35527
Calculated MW	62kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000FC:1/200 - 1/400ELISA:1/10000


Images



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Western blot analysis using KRT9 mouse mAb against HL-60 (1), A431 (2), Hela (3), and COS7 (4) cell lysate.



Flow cytometric analysis of HeLa cells using KRT9 mouse mAb (green) and negative control (red).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.