

KSHV K8a Mouse mAb

Catalog No: #63853

Package Size: #63853-1 50ul #63853-2 100ul


Description

Product Name	KSHV K8a Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG2a
Applications	WB;ICC
Immunogen Description	Purified recombinant fragment of KSHV K8α expressed in E. Coli.
Target Name	KSHV K8a
Other Names	K8
Formulation	Ascitic fluid containing 0.03% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

ICC:1/25-1/50ELISA:1/10000

Images



Western blot analysis using KSHV K8 mouse mAb against BCBL-1 (1) and TPA induced BCBL-1 (2) cell lysate.



Immunocytochemistry analysis of TPA induced BCBL-1 cells (A) and uninduced BCBL-1 cells (B) using KSHV K8 mouse mAb with AEC staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.