

KSHV ORF45 Mouse mAb

Catalog No: #63855

Package Size: #63855-1 50ul #63855-2 100ul


Description

Product Name	KSHV ORF45 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	ICC
Species Reactivity	KSHV
Immunogen Description	Purified recombinant fragment of KSHV ORF45 expressed in E. Coli.
Target Name	KSHV ORF45
Other Names	ORF45
Formulation	Purified antibody in PBS with 0.05% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

IHC:1/200 - 1/1000ELISA:1/10000

Images



Immunocytochemistry analysis of TPA induced BCBL-1 cells (A) and uninduced BCBL-1 cells (B) using KSHV ORF45 mouse mAb with AEC staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.