

LAMB1 Mouse mAb

Catalog No: #63857

Package Size: #63857-1 50ul #63857-2 100ul


Description

Product Name	LAMB1 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG2a
Applications	IHC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of LAMB1 (aa31-270) expressed in E. Coli.
Target Name	LAMB1
Other Names	CLM; MGC142015
Accession No.	P07942
Formulation	Purified antibody in PBS with 0.05% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

IHC:1/200 - 1/1000ELISA:1/10000

Images



Immunohistochemical analysis of paraffin-embedded human Thyroid tissues using LAMB1 mouse mAb

Note: This product is for in vitro research use only and is not intended for use in humans or animals.