

LCN1 Mouse mAb

Catalog No: #63859

Package Size: #63859-1 50ul #63859-2 100ul


Description

Product Name	LCN1 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	IHC;IF
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of LCN1 expressed in E. Coli.
Target Name	LCN1
Other Names	TP; PMFA; VEGP
Accession No.	P31025
Calculated MW	19.2kDa
Formulation	Ascitic fluid containing 0.03% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

IHC:1/200 - 1/1000ICC:1/200 - 1/1000ELISA:1/10000

Images



Confocal Immunofluorescence analysis of methanol-fixed HEK293 cells transfected with LCN1-hlgGfc using LCN1 mouse mAb(green), showing membrane localization. Blue: DRAQ5 fluorescent DNA dye.



Immunohistochemical analysis of paraffin-embedded human Prostate tissues using LCN1 mouse mAb

Note: This product is for in vitro research use only and is not intended for use in humans or animals.