

LDLR Mouse mAb

Catalog No: #63860

Package Size: #63860-1 50ul #63860-2 100ul


Description

Product Name	LDLR Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	IHC;FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human LDLR (AA: 22-150) expressed in E. Coli.
Target Name	LDLR
Other Names	FH; FHC; LDLCQ2
Accession No.	P01130
Calculated MW	95.4kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

IHC:1/200 - 1/1000FC:1/200 - 1/400ELISA:1/10000


Images



Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);



Flow cytometric analysis of Hela cells using LDLR mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human ovarian cancer tissues using LDLR mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.