

LEF1 Mouse mAb

Catalog No: #63861

Package Size: #63861-1 50ul #63861-2 100ul


Description

Product Name	LEF1 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;IHC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human LEF1 (AA: 33-138) expressed in E. Coli.
Target Name	LEF1
Other Names	LEF-1; TCF10; TCF7L3; TCF1ALPHA
Accession No.	Q9UJU2
Calculated MW	44.2kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000IHC:1/200 - 1/1000ELISA:1/10000


Images




Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);



Western blot analysis using LEF1 mouse mAb against Jurkat (1), HepG2 (2), and MOLT4 (3) cell lysate.



Immunohistochemical analysis of paraffin-embedded human colon cancer tissues using LEF1 mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human esophageal cancer tissues using LEF1 mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.