

LGR5 Mouse mAb

Catalog No: #63863

Package Size: #63863-1 50ul #63863-2 100ul


Description

Product Name	LGR5 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG2b
Applications	FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human LGR5 (AA: 22-178) expressed in E. Coli.
Target Name	LGR5
Other Names	FEX; HG38; GPR49; GPR67; GRP49
Accession No.	O75473
Calculated MW	100kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

FC:1/200 - 1/400ELISA:1/10000

Images



Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);



Flow cytometric analysis of A549 cells using LGR5 mouse mAb (green) and negative control (red).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.