

## LILRA3 Mouse mAb

Catalog No: #63867

Package Size: #63867-1 50ul #63867-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

|                       |                                                                                          |
|-----------------------|------------------------------------------------------------------------------------------|
| Product Name          | LILRA3 Mouse mAb                                                                         |
| Host Species          | Mouse                                                                                    |
| Clonality             | Monoclonal                                                                               |
| Isotype               | Mouse IgG1                                                                               |
| Applications          | FC                                                                                       |
| Species Reactivity    | Human                                                                                    |
| Immunogen Description | Purified recombinant fragment of human LILRA3 (AA: 24-168) expressed in E. Coli.         |
| Target Name           | LILRA3                                                                                   |
| Other Names           | HM31; HM43; ILT6; LIR4; CD85E; ILT-6; LIR-4                                              |
| Accession No.         | Q8N6C8                                                                                   |
| Calculated MW         | 47.5kDa                                                                                  |
| Formulation           | Purified antibody in PBS with 0.05% sodium azide                                         |
| Storage               | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles. |


## Application Details

FC:1/200 - 1/400ELISA:1/10000

## Images



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Flow cytometric analysis of K562 cells using LILRA3 mouse mAb (green) and negative control (red).

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.