

LILRB2 Mouse mAb

Catalog No: #63871

Package Size: #63871-1 50ul #63871-2 100ul


Description

Product Name	LILRB2 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG2a
Applications	WB;FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human LILRB2 (AA: 51-184) expressed in E. Coli.
Target Name	LILRB2
Other Names	ILT4; LIR2; CD85D; ILT-4; LIR-2; MIR10; MIR-10
Accession No.	Q8N423
Calculated MW	65kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000FC:1/200 - 1/400ELISA:1/10000


Images



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Western blot analysis using LILRB2 mouse mAb against MOLT4 (1) and K562 (2) cell lysate.



Flow cytometric analysis of HL-60 cells using LILRB2 mouse mAb (green) and negative control (red).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.