

LMO2 Mouse mAb

Catalog No: #63875

Package Size: #63875-1 50ul #63875-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	LMO2 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human LMO2 (AA: 1-158) expressed in E. Coli.
Target Name	LMO2
Other Names	TTG2; RBTN2; RHOM2; RBTNL1
Accession No.	P25791
Calculated MW	18.4kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

FC:1/200 - 1/400ELISA:1/10000

Images



Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);

Note: This product is for in vitro research use only and is not intended for use in humans or animals.