

LPL Mouse mAb

Catalog No: #63877

Package Size: #63877-1 50ul #63877-2 100ul


Description

Product Name	LPL Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of LPL expressed in E. Coli.
Target Name	LPL
Other Names	LIPD, HDLCQ11
Accession No.	P06858
Calculated MW	53.1kDa
Formulation	Ascitic fluid containing 0.03% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

WB:1/500 - 1/2000ELISA:1/10000

Images



Western blot analysis using LPL mouse mAb against Hela cell lysate (1).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.