

LRP1 Mouse mAb

Catalog No: #63879

Package Size: #63879-1 50ul #63879-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	LRP1 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	IHC;IF;FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human LRP1 (AA: 20-155) expressed in E. Coli.
Target Name	LRP1
Other Names	APR; LRP; A2MR; CD91; APOER; LRP1A; TGFBR5; IGFBP3R
Accession No.	Q07954
Calculated MW	50.5kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

IHC:1/200 - 1/1000 ICC:1/200 - 1/1000 FC:1/200 - 1/400 ELISA:1/10000


Images




Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);



Flow cytometric analysis of Jurkat cells using LRP1 mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human esophageal cancer tissues using LRP1 mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human stomach cancer tissues using LRP1 mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.