

LRP4 Mouse mAb

Catalog No: #63883

Package Size: #63883-1 50ul #63883-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	LRP4 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG2a
Applications	FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human LRP4 (AA: extra 21-189) expressed in E. Coli.
Target Name	LRP4
Other Names	CLSS; CMS17; LRP-4; LRP10; MEGF7; SOST2
Accession No.	O75096
Calculated MW	212kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

FC:1/200 - 1/400ELISA:1/10000

Images



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Flow cytometric analysis of HL-60 cells using LRP4 mouse mAb (green) and negative control (red).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.