

LRP5 Mouse mAb

Catalog No: #63884

Package Size: #63884-1 50ul #63884-2 100ul


Description

Product Name	LRP5 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human LRP5 expressed in E. Coli.
Target Name	LRP5
Other Names	HBM; LR3; OPS; EVR1; EVR4; LRP7; OPPG; BMND1; OPTA1; VBCH2
Accession No.	O75197
Calculated MW	179kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

FC:1/200 - 1/400ELISA:1/10000

Images



Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);

Note: This product is for in vitro research use only and is not intended for use in humans or animals.