

LSD1/AOF2 Mouse mAb

Catalog No: #63889

Package Size: #63889-1 50ul #63889-2 100ul


Description

|                       |                                                                                          |
|-----------------------|------------------------------------------------------------------------------------------|
| Product Name          | LSD1/AOF2 Mouse mAb                                                                      |
| Host Species          | Mouse                                                                                    |
| Clonality             | Monoclonal                                                                               |
| Isotype               | Mouse IgG1                                                                               |
| Applications          | WB;IHC                                                                                   |
| Species Reactivity    | Human,Mouse,Monkey                                                                       |
| Immunogen Description | Purified recombinant fragment of human LSD1 expressed in E. Coli.                        |
| Target Name           | LSD1/AOF2                                                                                |
| Other Names           | KDM1; LSD1; AOF2                                                                         |
| Accession No.         | O60341                                                                                   |
| Calculated MW         | 93kDa                                                                                    |
| Formulation           | Purified antibody in PBS with 0.05% sodium azide.                                        |
| Storage               | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles. |


Application Details

WB:1/500 - 1/2000IHC:1/200 - 1/1000ELISA:1/10000

Images



Western blot analysis using LSD1 mouse mAb against COS (1), Hela (2), NIH/3T3 (3), A549 (4) and Jurkat (5) cell lysate.



Immunohistochemical analysis of paraffin-embedded human lung carcinoma (left) and kidney carcinoma (right), showing nuclear localization using LSD1 mouse mAb with DAB staining.

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.