

LTBR Mouse mAb

Catalog No: #63890

Package Size: #63890-1 50ul #63890-2 100ul


Description

Product Name	LTBR Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG2a
Applications	WB;FC
Species Reactivity	Human,Mouse
Immunogen Description	Purified recombinant fragment of human LTBR (AA: extra 31-227) expressed in E. Coli.
Target Name	LTBR
Other Names	TNFCR; TNFR3; D12S370; TNFR-RP; TNFRSF3; TNFR2-RP; LT-BETA-R; TNF-R-III
Accession No.	P36941
Calculated MW	46.7kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details


WB:1/500 - 1/2000FC:1/200 - 1/400ELISA:1/10000

Images



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)





Flow cytometric analysis of HL-60 cells using LTBR mouse mAb (green) and negative control (red).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.