

LYN Mouse mAb

Catalog No: #63891

Package Size: #63891-1 50ul #63891-2 100ul


Description

Product Name	LYN Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG2b
Applications	WB
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of LYN expressed in E. Coli.
Target Name	LYN
Other Names	JTK8; FLJ26625
Accession No.	P07948
Calculated MW	56kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

WB:1/500 - 1/2000ELISA:1/10000

Images



Western blot analysis using LYN mouse mAb againsts HL60 (1), L540 (2), SLLP-M2 (3), SEM (4) and Ramos (5) cell lysate.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.