

M6PR Mouse mAb

Catalog No: #63892

Package Size: #63892-1 50ul #63892-2 100ul


Description

Product Name	M6PR Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG2b
Applications	WB;IHC;FC
Species Reactivity	Human,Mouse
Immunogen Description	Purified recombinant fragment of human M6PR (AA: 124-277) expressed in E. Coli.
Target Name	M6PR
Other Names	SMPR; MPR46; CD-MPR; MPR 46; MPR-46; CD-M6PR
Accession No.	P20645
Calculated MW	30.9kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000IHC:1/200 - 1/1000FC:1/200 - 1/400ELISA:1/10000


Images




Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)




Western blot analysis using M6PR mouse mAb against mouse brain (1) and HepG2 (2) cell lysate.



Flow cytometric analysis of THP-1 cells using M6PR mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human ovarian cancer tissues using M6PR mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human brain tissues using M6PR mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.